



Product Manual

OxiSelect™ N^ε-(carboxymethyl) lysine (CML) Competitive ELISA Kit

Catalog Number

STA-816	96 assays
STA-816-5	5 x 96 assays

**NOTE: Revisions to
"Preparation of Reagents"**

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures

Gentaur Molecular Products BVBA
Address: Voortstraat 49, 1910 Kampenhout, Belgium
T: 0032 16 58 90 45 | **E:** info@gentaur.com
Websites: www.gentaur.com | www.maxanim.com

Introduction

The non-enzymatic reaction of reducing carbohydrates with lysine side chains and N-terminal amino groups of macromolecules (proteins, phospholipids and nucleic acids) is called the Maillard reaction or glycation. The products of this process, termed advanced glycation end products (AGEs), adversely affect the functional properties of proteins, lipids and DNA. Tissue levels of AGE increase with age and the formation of AGEs is predominantly endogenous, though these products can also be derived from exogenous sources such as food and tobacco smoke. AGE modification of proteins can contribute to the pathophysiology of aging and long-term complications of diabetes, atherosclerosis and renal failure. AGEs also interact with a variety of cell-surface AGE-binding receptors (RAGE), leading either to their endocytosis and degradation or to cellular activation and pro-oxidant or pro-inflammatory events.

Although several AGE structures have been reported, it was demonstrated that N^ε-(carboxymethyl) lysine (CML) is a major antigenic AGE structure. CML concentration is increased in patients who have diabetes with complications, including nephropathy, retinopathy, and atherosclerosis. CML is also recognized by receptor for AGE (RAGE), and CML-RAGE interaction activates cell signaling pathways such as NF-κB.




Figure 1. Structure of N^ε-(carboxymethyl) lysine (CML)

OxiSelect™ N^ε-(carboxymethyl) lysine (CML) Competitive ELISA Kit provides rapid detection and quantitation of CML protein adducts. The quantity of CML adduct in protein samples is determined by comparing its absorbance with that of a known CML-BSA standard curve. Each kit provides sufficient reagents to perform up to 96 assays, including standard curve and unknown protein samples.

Assay Principle

First, a CML conjugate is coated on the ELISA plate. The unknown CML protein samples or CML-BSA standards are then added to the CML conjugate preabsorbed plate. After a brief incubation, the anti-CML monoclonal antibody is added, followed by an HRP conjugated secondary antibody. The content of CML protein adducts in unknown samples is determined by comparison with the predetermined CML-BSA standard curve.

Related Products

1. STA-305: OxiSelect™ Nitrotyrosine ELISA Kit
2. STA-310: OxiSelect™ Protein Carbonyl ELISA Kit

3. STA-811: OxiSelect™ Methylglyoxal (MG) Competitive ELISA Kit
4. STA-813: OxiSelect™ N^ε-(carboxyethyl) lysine (CEL) Competitive ELISA Kit
5. STA-817: OxiSelect™ Advanced Glycation End Products (AGE) Competitive ELISA Kit
6. STA-832: OxiSelect™ MDA Adduct Competitive ELISA Kit
7. STA-838: OxiSelect™ HNE Adduct Competitive ELISA Kit

Kit Components

Box 1 (shipped at room temperature)

1. 96-well Protein Binding Plate (Part No. 231001): One strip well 96-well plate.
2. Anti-CML Antibody (1000X) (Part No. 281601): One 10 µL vial of anti-CML antibody.
3. Secondary Antibody, HRP Conjugate (1000X) (Part No. 230003): One 20 µL vial.
4. Assay Diluent (Part No. 310804): One 50 mL bottle.
5. 10X Wash Buffer (Part No. 310806): One 100 mL bottle.
6. Substrate Solution (Part No. 310807): One 12 mL amber bottle.
7. Stop Solution (Part. No. 310808): One 12 mL bottle.

Box 2 (shipped on blue ice packs)

1. CML-BSA Standard (Part No. 231602): One 50 µL vial of 1.0 mg/mL CML-BSA in PBS. CML-BSA is prepared as described by Koito *et al.* (see Ref. 9) and it has 15 moles of CML per mole of BSA.
2. 1000X CML Conjugate (Part No. 281602): One 20 µL vial.
3. 100X Conjugate Diluent (Part No. 281603): One 300 µL vial.

Materials Not Supplied

1. Protein samples such as purified protein, plasma, serum, cell lysate
2. 1X PBS
3. 10 µL to 1000 µL adjustable single channel micropipettes with disposable tips
4. 50 µL to 300 µL adjustable multichannel micropipette with disposable tips
5. Multichannel micropipette reservoir
6. Microplate reader capable of reading at 450 nm (620 nm as optional reference wave length)

Storage

Aliquot and store the Anti-CML Antibody, CML-BSA Standard, CML Conjugate and 100X Conjugate Diluent at -20°C to avoid multiple freeze/thaw cycles. Store all other components at 4°C.

Preparation of Reagents

- CML Conjugate Coated Plate:

Note: The CML Conjugate coated wells are not stable and should be used within 24 hrs after coating. Only coat the number of wells to be used immediately.

1. Immediately before use, prepare 1X Conjugate Diluent by diluting the 100X Conjugate Diluent in 1X PBS. Example: Add 50 μ L to 4.95 mL of 1X PBS.
 2. Immediately before use, prepare 1X CML Conjugate by diluting the 1000X CML Conjugate in 1X Conjugate Diluent. Example: Add 5 μ L of 1000X CML Conjugate to 4.995 mL of 1X Conjugate Diluent.
 3. Add 100 μ L of the 1X CML Conjugate to each well to be tested and incubate overnight at 4°C. Remove the CML Conjugate coating solution and wash twice with 1X PBS. Blot plate on paper towels to remove excess fluid. Add 200 μ L of Assay Diluent to each well and block for 1 hr at room temperature on an orbital shaker. Transfer the plate to 4°C and remove the Assay Diluent **immediately before use**.
- 1X Wash Buffer: Dilute the 10X Wash Buffer Concentrate to 1X with deionized water. Stir to homogeneity.
 - Anti-CML Antibody and Secondary Antibody: Immediately before use, dilute the Anti-CML antibody 1:1000 and Secondary Antibody 1:1000 with Assay Diluent. Do not store diluted solutions.

Preparation of Standard Curve

Prepare a dilution series of CML-BSA standards in the concentration range of 0 to 12.5 μ g/mL by diluting the CML-BSA Standard in Assay Diluent (Table 1).

Standard Tubes	1 mg/mL CML-BSA Standard (μ L)	Assay Diluent (μ L)	CML-BSA (μ g/mL)	CML (ng/mL)
1	5	395	12.5	576
2	200 of Tube #1	200	6.25	288
3	200 of Tube #2	200	3.13	144
4	200 of Tube #3	200	1.56	72
5	200 of Tube #4	200	0.78	36
6	200 of Tube #5	200	0.39	18
7	200 of Tube #6	200	0.20	9
8	200 of Tube #7	200	0.10	4.5
9	200 of Tube #8	200	0.050	2.25
10	0	200	0	0

Table 1. Preparation of CML-BSA Standards

Assay Protocol

Note: If testing mouse or rat plasma or serum, the IgG must be completely removed from each sample prior to testing, such as with Protein A or G beads. Additionally, a control well without primary antibody should be run for each sample to determine background signal.

1. Prepare and mix all reagents thoroughly before use. Each CML sample including unknown and standard should be assayed in duplicate.
2. Add 50 μ L of unknown sample or CML-BSA standard to the wells of the CML Conjugate coated plate. If needed, unknown samples may be diluted in 1X PBS containing 0.1% BSA before adding. Incubate at room temperature for 10 minutes on an orbital shaker.
3. Add 50 μ L of the diluted anti-CML antibody to each well, incubate at room temperature for 1 hour on an orbital shaker.
4. Wash 3 times with 250 μ L of 1X Wash Buffer with thorough aspiration between each wash. After the last wash, empty wells and tap microwell strips on absorbent pad or paper towel to remove excess 1X Wash Buffer.
5. Add 100 μ L of the diluted Secondary Antibody-HRP Conjugate to all wells and incubate for 1 hour at room temperature on an orbital shaker. Wash the strip wells 3 times according to step 4 above.
6. Warm Substrate Solution to room temperature. Add 100 μ L of Substrate Solution to each well. Incubate at room temperature for 2-20 minutes on an orbital shaker.
Note: Watch plate carefully; if color changes rapidly, the reaction may need to be stopped sooner to prevent saturation.
7. Stop the enzyme reaction by adding 100 μ L of Stop Solution to each well. Results should be read immediately (color will fade over time).
8. Read absorbance of each well on a microplate reader using 450 nm as the primary wave length.

Example of Results

The following figures demonstrate typical CML Competitive ELISA results. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 2: CML-BSA Competitive ELISA Standard Curve.




Figure 3: CML Protein Adduct in Human Plasma.

References

1. Monnier, V., and Cerami, A. (1981) *Science* **211**, 491–493.
2. Ahmed M.U., Thorpe S.R., Baynes J.W (1986) *J. Biol. Chem.* **261**, 4889–4894.
3. Reddy S., Bichler J., Wells-Knecht K.J., Thorpe S.R., Baynes J.W (1995) *Biochemistry* **34**, 10872–10878.
4. Dunn, J. A., Patrick, J. S., Thorpe, S. R., and Baynes, J. W. (1989) *Biochemistry* **28**, 9464-9468.
5. Ahmed, M. U., Brinkmann Frye, E., Degenhardt, T. P., Thorpe, S. R., and Baynes, J. W. (1997) *Biochem. J.* **324**, 565-570.
6. Sell, D. R., and Monnier, V. M. (1989) *J. Biol. Chem.* **264**, 21597-21602.
7. Onorato, J., Jenkins, A., Thorpe, S., and Baynes, J. (2000) *J. Biol. Chem.* **275**, 21177–21184.
8. Boehm BO, Schilling S, Rosinger S, Lang GE, Lang GK, Kientsch-Engel P, Stahl P (2004) *Diabetologia* **47**, 1376–1379.
9. Koito W, Araki T, Horiuchi S, Nagai R (2004) *J. Biochem.* **136**, 831-837.

Recent Product Citations

1. Shen, C.Y. et al. (2023). Unveiling the molecular basis of inflamm-aging induced by advanced glycation end products (AGEs)-modified human serum albumin (AGE-HSA) in patients with different immune-mediated diseases. *Clin Immunol.* **252**:109655. doi: 10.1016/j.clim.2023.109655.
2. Damasiewicz-Bodzek, A. & Nowak, A. (2022). Concentrations of N6-Carboxymethyllysine (CML), N6-Carboxyethyllysine (CEL), and Soluble Receptor for Advanced Glycation End-Products (sRAGE) Are Increased in Psoriatic Patients. *Biomolecules.* **12**(12):1870. doi: 10.3390/biom12121870.
3. Vaidya, R. et al. (2022). Accumulation of fluorescent advanced glycation end products and carboxymethyl-lysine in human cortical and trabecular bone. *Bone Rep.* doi: 10.1016/j.bonr.2022.101634.
4. Imai, Y. et al. (2022). Inhibitory Effects of Parachlorella Beijerinckii Extracts on the Formation of Advanced Glycation End Products and Glycative Stress-Induced Inflammation in an In Vitro Skin Dermis-Like Model. *Evid. Based Complementary Altern. Med.* doi: 10.1155/2022/8789903.
5. Yoon, S. et al. (2022). Effect of Cirsium japonicum Flower Extract on Skin Aging Induced by Glycation. *Molecules.* **27**(7):2093. doi: 10.3390/molecules27072093.
6. Maciejczyk, M. et al. (2022). Oxidation, Glycation, and Carbamylation of Salivary Biomolecules in Healthy Children, Adults, and the Elderly: Can Saliva Be Used in the Assessment of Aging? *J Inflamm Res.* **15**:2051-2073. doi: 10.2147/JIR.S356029.
7. Jarrete, A.P. et al. (2022). Alterations in pro- and anti-inflammatory mediators are involved in microvascular dysfunction in postmenopausal women with type 2 diabetes mellitus. *Braz J Med Biol Res.* **55**:e11821. doi: 10.1590/1414-431X2021e11821.
8. Rajamanickam, A. et al. (2021). Diminished circulating levels of angiogenic factors and RAGE ligands in helminth-diabetes comorbidity and reversal following anthelmintic treatment. *J Infect Dis.* doi: 10.1093/infdis/jiab170.
9. Pantner, Y. et al. (2021). DJ-1 attenuates the glycation of mitochondrial complex I and complex III in the post-ischemic heart. *Sci Rep.* **11**(1):19408. doi: 10.1038/s41598-021-98722-1.
10. Shin, S. et al. (2021). Anti-glycation activities of methyl gallate in vitro and in human explants. *J Cosmet Dermatol.* doi: 10.1111/jocd.14406.
11. Mahmoud, A.M. & Ali, M.M. (2021) High Glucose and Advanced Glycation End Products Induce CD147-Mediated MMP Activity in Human Adipocytes. *Cells.* **10**(8):2098. doi: 10.3390/cells10082098.
12. Lazzari, T.K. et al. (2021). Leptin and advanced glycation end products receptor (RAGE) in tuberculosis patients. *PLoS One.* **16**(7):e0254198. doi: 10.1371/journal.pone.0254198.
13. Ahmad, S. et al. (2021). Gold Nanoparticle-Bioconjugated Aminoguanidine Inhibits Glycation Reaction: An In Vivo Study in a Diabetic Animal Model. *Biomed Res Int.* doi: 10.1155/2021/5591851.
14. Chiew, Y. et al. (2021). Tocotrienol-rich vitamin E from palm oil (Tocovid) and its effects in diabetes and diabetic retinopathy: a pilot phase II clinical trial. *Asian J. Ophthalmol.* **17**(4):375-399. doi: 10.35119/asjoo.v17i4.698.
15. Altomare, A. et al. (2021). In-Depth AGE and ALE Profiling of Human Albumin in Heart Failure: Ex Vivo Studies. *Antioxidants (Basel).* **10**(3):358. doi: 10.3390/antiox10030358.
16. Li, Y.Y. et al. (2021). Protective effects of dietary carnosine during in-vitro digestion of pork differing in fat content and cooking conditions. *J Food Biochem.* **45**(2):e13624. doi: 10.1111/jfbc.13624.

17. Chen, Z. et al. (2020). Association of carbamylated high-density lipoprotein with coronary artery disease in type 2 diabetes mellitus: carbamylated high-density lipoprotein of patients promotes monocyte adhesion. *J Transl Med.* **18**(1):460. doi: 10.1186/s12967-020-02623-2.
18. Merhi, Z. et al. (2020). Perinatal Exposure to High Dietary Advanced Glycation End-Products Affects the Reproductive System in Female Offspring in Mice. *Mol Hum Reprod.* doi: 10.1093/molehr/gaaa046.
19. Gutierrez-Mariscal, F.M. et al. (2020). Reduction in Circulating Advanced Glycation End Products by Mediterranean Diet is Associated with Increased Likelihood of type 2 Diabetes Remission in Patients with Coronary Heart Disease: From the Cordioprev Study. *Mol Nutr Food Res.* doi: 10.1002/mnfr.201901290.
20. Thornton, K. et al. (2020). Dietary Advanced Glycation End Products (AGEs) could alter ovarian function in mice. *Mol Cell Endocrinol.* doi: 10.1016/j.mce.2020.110826.
21. Hernández, C. et al. (2020). The Usefulness of Serum Biomarkers in the Early Stages of Diabetic Retinopathy: Results of the EUROCONDOR Clinical Trial. *J Clin Med.* **9**(4). pii: E1233. doi: 10.3390/jcm9041233.
22. Velayoudom-Cephise, F.L. et al. (2020). Receptor For Advanced Glycated End Products Modulates Oxidative Stress And Mitochondrial Function In The Soleus Muscle Of High Fat Fed Mice. *Appl Physiol Nutr Metab.* doi: 10.1139/apnm-2019-0936.
23. Chen, S.H. et al. (2020). Iron and Advanced Glycation End Products: Emerging Role of Iron in Androgen Deficiency in Obesity. *Antioxidants.* **9**:261. doi: 10.3390/antiox9030261.
24. Shimizu, Y. et al. (2020). Role of DJ-1 in Modulating Glycative Stress in Heart Failure. *J Am Heart Assoc.* **9**(4). doi: 10.1161/jaha.119.014691.
25. de la Cruz-Ares, S. et al. (2020). Endothelial Dysfunction and Advanced Glycation End Products in Patients with Newly Diagnosed Versus Established Diabetes: From the CORDIOPREV Study. *Nutrients.* **12**(1). pii: E238. doi: 10.3390/nu12010238.
26. Lee, J. et al. (2019). Mitochondrial carnitine palmitoyltransferase 2 is involved in N ϵ -(carboxymethyl)-lysine-mediated diabetic nephropathy. *Pharmacol Res.* doi: 10.1016/j.phrs.2019.104600.
27. Kaburagi, T. et al. (2019). Low-Carbohydrate Diet Inhibits Different Advanced Glycation End Products in Kidney Depending on Lipid Composition but Causes Adverse Morphological Changes in a Non-Obese Model Mice. *Nutrients.* **11**(11). pii: E2801. doi: 10.3390/nu11112801.
28. Yang, J. et al. (2019). Neutrophil-derived advanced glycation end products-N ϵ -(carboxymethyl) lysine promotes RIP3-mediated myocardial necroptosis via RAGE and exacerbates myocardial ischemia/reperfusion injury. *FASEB J.* doi: 10.1096/fj.201900115RR.
29. Ndidi, U.S. et al. (2019). Effect of N(Epsilon)-(carboxymethyl)lysine on Laboratory Parameters and Its Association with β S Haplotype in Children with Sickle Cell Anemia. *Disease Markers.* doi: 10.1155/2019/1580485.
30. Ferron, A.J.T. et al. (2019). Protective Effect of Tomato-Oleoresin Supplementation on Oxidative Injury Recoveries Cardiac Function by Improving β -Adrenergic Response in a Diet-Obesity Induced Model. *Antioxidants (Basel).* **8**(9). pii: E368. doi: 10.3390/antiox8090368.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products.

In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.